

A photograph of a person lying on their back on a dark asphalt surface. The person is wearing a blue and black plaid shirt, tan trousers, blue socks, and black sneakers. Their arms are outstretched to the sides, and their legs are bent at the knees. A banana peel is visible on the ground near their feet. The background is a solid blue color with a white geometric pattern on the left side.

# Good bye slippage

– a new fusion to tackle bandage slippage on the foot

Presenters:

Josefin Damm & Andreas Nilsson



# Disclaimer/Conflict of Interest

- Josefin Damm: Co-inventor of Lundatex<sup>®</sup> products, co-founder of PressCise
- Andreas Nilsson: CEO of PressCise
- None of the presenters are nurse or MD

# Background

“In fact, the medical field of compression treatment is maybe the only one where quantitative dosage has almost never been measured, despite outcomes largely depending on it.”



Mosti & Partsch,  
Eur J Vasc Endovasc Surg (2017)



# Problem with bandage application – no consistency in pressure

Protz et.al. 2014

- 551 healthcare personnel
- Target pressure: 50-60 mmHg





# Our solution – to control the pressure



Precise pressure over the entire leg:

- Invariant of different appliers
- Different sizes
- Different shapes

(Oedema – swelling/deswelling)


# Study by Wiklander et.al. (2015)

International Wound Journal ISSN 1742-4801

## ORIGINAL ARTICLE

### An investigation of the ability to produce a defined 'target pressure' using the PressCise compression bandage

Kerstin Wiklander<sup>1</sup>, Annette Erichsen Andersson<sup>2</sup> & Ulrika Källman<sup>3,4</sup>



# Study by G.Mosti &H.Partsch (2017)

Eur J Vasc Endovasc Surg (2017) ■, 1–6

## **A New Two Component Compression System Turning an Elastic Bandage into an Inelastic Compression Device: Interface Pressure, Stiffness, and Haemodynamic Effectiveness**

Giovanni Mosti <sup>a,\*</sup>, Hugo Partsch <sup>b</sup>

<sup>a</sup>Department of Angiology, MD Barbantini Clinic, Lucca, Italy

<sup>b</sup>Professor Emeritus Medical University of Vienna, Vienna, Austria

- 25 legs from 25 patients
- Venous insufficiency

*“...all affected by clinically significant reflux in the great saphenous vein (GSV), with clinical stage C2-C5”*



# Study by G.Mosti &H.Partsch (2017)

Eur J Vasc Endovasc Surg (2017) ■, 1–6

**A New Two Component Compression System Turning an Elastic Bandage into an Inelastic Compression Device: Interface Pressure, Stiffness, and Haemodynamic Effectiveness**

Giovanni Mosti <sup>a,\*</sup>, Hugo Partsch <sup>b</sup>

<sup>a</sup>Department of Angiology, MD Barbantini Clinic, Lucca, Italy  
<sup>b</sup>Professor Emeritus Medical University of Vienna, Vienna, Austria




Interface pressure

B1-point: Medial side gastrocnemius muscle turns into the tendinous part


C-point: at maximum calf circumference

supine and standing pressure  
Lundatex 20



	supine	standing
Minimum	16.00	17.00
Median	20.00	21.00
Maximum	25.00	32.00
Mean	20.40	22.00
Std. Deviation	2.179	3.674

supine and standing pressure  
Lundatex 20



	supine	standing
Minimum	16.00	17.00
Median	23.00	25.00
Maximum	28.00	37.00
Mean	23.28	24.80
Std. Deviation	2.670	4.320


# Uniform & constant pressure



# How does it work?

- Lundatex medical is based on Laplace's law

Pressure = force × overlap × curvature







# How does it work?



- Lundatex medical is based on Laplace's law

Pressure = force × **overlap** × curvature







$$\text{Pressure} = \text{force} \times \text{overlap} \times \text{curvature}$$



# Bandages with markings

- same same but different?



$$\text{Pressure} = \text{force} \times \text{overlap} \times \text{curvature}$$





# Why we need innovation in stockings

- Fitting / Custom made stockings
  - Pressure?
- Oedema – swelling/deswelling
  - Pressure?
- Comfort
  - patient compliance
- Donning problems
  - varies with compression class and elasticity of the material used



Reality

# A new Smart Textile Stocking



New properties allows for:

- Well-defined pressure (20 mmHg)
- Uniform pressure
- Ensured pressure regardless of leg shape
- Easier donning

# Two limited pilot studies – methods

1. 1<sup>st</sup> study: 10 healthy subjects (5 women and 3 men)
2. 2<sup>nd</sup> study: 8 healthy subjects (5 women and 5 men)

Interface pressure was measured with a Picopress

At B1-level gastrocnemius muscle turns into the tendinous part:

B1 – medial

B2 – lateral

At C level - point at maximum leg-circumference:

C1 – medial

C2 – posterior


C3 – lateral





# Results – 1<sup>st</sup> pilot study

- Leg circumference
- range 20.7 – 29.3 cm



- range 30.3 – 51.0 cm

*Lundatex*<sup>®</sup> stocking


andreas@presscise.com

10 subjects - Supine		
Sensor Position	Interface pressure Mean (mmHg)	SD
B1	20.7	1.6
B2	21.1	1.7
C1	20.9	1.7
C2	21.9	1.7
C3	21.8	2.2

josefin@presscise.com

# Results – 2<sup>nd</sup> pilot study

*Same stocking used on all 8 subjects!*




*Lundatex*<sup>®</sup> stocking

andreas@presscise.com

josefin@presscise.com

# Uniform & constant pressure



# Problem with bandage application on the foot



Source: [www.worlwidewounds.com](http://www.worlwidewounds.com)

*“Bandage slippage can create local high-pressure areas that may cause tissue damage and even necrosis.”*

S. Thomas, “Compression Bandaging in the Treatment of Venous Leg Ulcers,” *World Wide Wounds*, 1998.



# What if we take the best from two worlds and combine into one solution?


Potential dangerous zone!

- Too low pressure
- Too high pressure



# Pilot study, method – sock-bandage fusion

- Thirteen healthy subjects (six females, seven males)
- Three pressure sensors were placed unilaterally:  
(A) on the foot, (B) at the ankle and (C) on the calf



# Pilot study, method – sock-bandage fusion

1. Sock (20 mmHg) covering sensor A
2. 5 cm cuff of the sock (10 mmHg) covering sensor B
3. 5 cm first bandage turn (10 mmHg) covering sensor B
4. Bandage (20 mmHg) covering sensor C
5. The interface pressure was measured with a Picopress® in supine and standing



# Results – sock-bandage fusion



*Same sock used on all 13 subjects!*

	Foot (A)		Ankle (B)		Leg (C)	
	Mean	SD	Mean	SD	Mean	SD
Supine	<b>20</b>	2.9	<b>20</b>	2.4	<b>21</b>	1.6
Standing	<b>20</b>	3.9	<b>21</b>	3.2	<b>22</b>	2.0
	Sock		Transition area		Bandage	

*The sock-bandage fusion method applied a well-defined pressure at the foot, ankle and leg*



# Results – sock-bandage fusion



# Conclusion – sock-bandage fusion

- Well defined pressure – on foot, ankle and leg
- Easy donning – on the foot and leg
- Comfortable – even pressure, no slippage
- Preserves normal ankle range-of-motion
- Ability to wear normal footwear – no bulky material



Good bye slippage!

